



# BIORAD®

## THERMO CLEAN PLUS

Hemodialysis Device  
Disinfection



### FORMS OF PACKAGE

Description	Package
5 L PE drum	4x5 L



### CHEMICAL COMPOSITION:

Acetic acid<%12 , peracetic acid<%20 , hydrogen peroxide< %20 and auxiliary agents

### PRODUCT SPECIFICATIONS:

THERMO CLEAN Plus provides cleaning, decalcification and disinfection of organic and inorganic waste that forms on the generator and lines of haemodialysis equipment during haemodialysis process. It is compliant with stainless steel, plastic materials and ceramic surfaces. It should not contact with plexiglass and aluminium surfaces

### INTENDED USE:

Thermo Clean Plus is compliant with stainless steel, glass, plastic materials and ceramic surfaces. Do not use it on sensitive surfaces such as plexiglass and aluminium that may react with acids.

### USAGE:

Thermo Clean Plus is used by automatic vacuum through disinfection pipette of the haemodialysis machines. Observe recommendations of the machine manufacturer. pH=1-3

### ANTIMICROBIAL CHARACTERISTICS:

Bactericidal, Fungicidal, Virucidal, Micobactericidal,Tuberculocidal