



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.

Features












- 10 x Magnification
- 32mm wide field of view
- Polarized, non-polarized, and amber light illumination
- Detachable protective glass
- Automatic shutdown
- Included adapter fits all phone
- All-metal housing

DERMATOSKOP PRO DE-3100



DERMATOSKOP PRO DE-3100

Technical Information

	Material	Aluminum & Glass
	Optical Design	All glass, 4 elements 3 groups, anti-reflection coating
	Lens Diameter	32mm(front); 25mm(rear)
	Magnification	10X
	Distortion	8%
	Polarization	Polarization & non-polarization & UV 4 pcs UV LED 15 pcs large size
	Resolution	300 LP/MM (Axis) 250 LP/MM (Edge)
	LED Type	SMD LED beads
	Battery Capacity	1000mAh Lithium ion
	Charging	USB-C
	Working Time	3-6 hours



www.famed.com.tr



info@famed.com.tr



Accessories

DE-3100



UNIVERSAL PHONE ADAPTER

Save images of lesions on your smartphone



IN THE BOX

Relative accessories included in the package



SMALL CONTACT PLATE

To detect difficulty access pigmented lesions




CAMERA ADAPTER

Save images of lesions on your camera



PHONE COVER ADAPTER

iPhone & Samsung phone cover for dermatoscope



BACK UP BATTERY

1000mAh Lithium ion battery for replacement

