



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.

DERMATOSCOPE DE-200



Features

- 6 x Magnification
- 20mm wide field of view
- Polarized light illumination
- Automatic shutdown
- Universal Phone Adapter
- All-metal housing



*Like a scalpel, the dermatoscope is an awesome instrument in trained hands.
It's also a great tool for self-diagnostics with the help of dermatologist or dermascope App.*

Technical Information



Material	Aluminum & Glass
Optical Design	All glass, 2 elements 2 groups, anti-reflection coating
Lens Diameter	20mm
Magnification	6X
Distortion	8%
Polarization	2 LEDs white polarized
Resolution	300 LP/MM (Axis) 250 LP/MM (Edge)
LED Type	SMD LED beads
Battery Capacity	200mAh Lithium ion
Charging	USB-C
Working Time	90 minutes

Address:: Cevizli Mh. Hacilar Cd. No:44/A Kartal - İstanbul / TURKEY
Phone : +90 216 597 0 567 / +90 216 441 29 99 Fax: +90 216 397 15 15



www.famed.com.tr



info@famed.com.tr