



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.

DERMATOSCOPE DE-300

Features

- 6 x Magnification
- 20mm wide field of view (Aperture)
- Polarized light illumination
- Automatic shutdown
- Universal Phone Adapter
- All-metal housing



Like a scalpel, a dermatoscope is a fantastic tool in trained hands. It's also a great tool for self-diagnosis with the help of a dermatologist or dermascope application.

Technical Information

	Material	Aluminum & Glass
	Optical Design	All glass, 2 elements 2 groups, anti-reflection coating
	Lens Diameter	30mm
	Magnification	6X
	Distortion	8%
	Polarization	2 LEDs white non polarized, 4 LEDs white polarized
	Resolution	300 LP/MM (Axis) 250 LP/MM (Edge)
	LED Type	SMD LED beads
	Battery Capacity	200mAh Lithium ion
	Charging	USB-C
	Working Time	90 minutes

Address:: Cevizli Mh. Hacilar Cd. No:44/A Kartal - İstanbul / TURKEY
Phone : +90 216 597 0 567 / +90 216 441 29 99 **Fax:** +90 216 397 15 15

www.famed.com.tr

info@famed.com.tr