



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.


Özellikleri

- 10 x Büyütme
- 32 mm geniş görüş alanı
- Polarize, polarize olmayan ve amber ışık aydınlatma
- Çıkarılabilir koruyucu cam
- Otomatik kapanma
- Birlikte verilen adaptör tüm telefonlara uyar
- Tamamı metal gövde

DERMATOSKOP DE-3100



DERMATOSKOP DE-3100



Teknik Bilgiler

Malzeme	Alüminyum ve Cam
Optik Tasarım	Tamamı cam, 4 elemanlı 3 grup, yansıma önleyici kaplama
Mercek Çapı	32 mm (ön); 25mm(arka)
Büyütme	10X
Çarpıtma	8%
Polarizasyon	4 LED beyaz polarizesiz, 4 LED sarı polarize, 6 adet beyaz sarı. 18 adet beyaz. 6 adet turuncu toplam 30 led
Çözünürlük	300 LP/MM (Eksen) 250 LP/MM (Kenar)
Led Tipi	SMD LED boncuklar
Pil Kapasitesi	1000mAh Lityum iyon
Dolum	USB-C
Çalışma Zamanı	3-6 Saat



www.famed.com.tr



info@famed.com.tr



Accessories

DE-3100



UNIVERSAL PHONE ADAPTER

Save images of lesions on your smartphone



IN THE BOX

Relative accessories included in the package



SMALL CONTACT PLATE

To detect difficulty access pigmented lesions



CAMERA ADDAPTER

Save images of lesions on your camera



PHONE COVER ADAPTER

iPhone & Samsung phone cover for dermatoscope



BACK UP BATTERY

1000mAh Lithium ion battery for replacement



Address:: Cevizli Mh. Hacilar Cd. No:44/A Kartal - İstanbul / TURKEY
Phone : +90 216 597 0 567 / +90 216 441 29 99 **Fax:** +90 216 397 15 15
W: www.famed.com.tr **E:** info@famed.com.tr