

AUTOMATIC ENDOSCOPE WASHER & DISINFECTOR
New Technology for Cleaning and Disinfection

FAMED SAĞLIK TIBBİ CİHAZLAR İTH.İHRC. ve SAN. TİC. LTD. ŞTİ.



*Automatic Endoscope Washer
Disinfector Device
Confidence from Medical Tests
Warranty.*

Choyang Medical Endoscope Washer-Disinfector offers all features necessary to provide triple protection - protecting patients, protecting staff, and protecting all kinds of your flexible endoscopes.

This product is certified as a medical device in the European Union under the Medical Device Directive 93/42/EEC by SGS CE0120, exclusively for the intended use of disinfects the flexible endoscope. Other non medical uses ascribed to this device are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those.

Multi-Channel jetting system rotated in a 360° arc



3 Adet güçlü motorla müttaşem yıkama ve dezenfeksiyon



Automatic Water Line Sterilization (self disinfection) cycle

UV lamp disinfection system in the water tank



SPECIFICATION

Leakage Tester (water Type) / Operating Procedure

Time	Washing	10~50sec, 1min~5min
	Disinfection	1 ~ 9 min, 10 ~ 50 min, 1 ~ 9 hr
	Enzimatik	10 sec
	Alcohol	5 sec
Program (Automatic)	Washing-Detergent-Washing-Disinfectant-Washing-Alcohol-Air	
TANK CAPACITY	Water	15 Liters
	Disinfectant	15 Liters
	Detergent :	1,5 Liters
	Alcohol :	1,5 Liters
COUPLERS	4 Outlets in cleaning Base (1 Scope)	
Discharge	1 Outlet (Compulsory discharging on rear panel)	
	3 Outlets (in the front door)	
Power	AC 220/230V, 50/60Hz	
Weight	550 x 560 x 940mm (W x D x H)	
Dimension	80 kg.	

FEATURES

- Single scope cleaning and disinfection
- Leakage tester
- Perfect washing and disinfection by 3 powerful motors
- Multi-channel jetting system rotated in a 360 arc
- Advanced water filtration system
- UV lamp water disinfection system
- Thermal control system for disinfection
- Self disinfection system
- Alcohol purging system
- Printer Connection (optional)