

AUTOMATIC ENDOSCOPE WASHER & DISINFECTOR
New Technology for Cleaning and Disinfection

FAMED SAĞLIK TIBBİ CİHAZLAR İTH.İHRC. ve SAN. TİC. LTD. ŞTİ.



**Automatic Endoscope Washer
Disinfecter Device**
Confidence from Medical Tests
Warranty.

Choyang Medical Endoscope Washer-Disinfecter offers all features necessary to provide triple protection - protecting patients, protecting staff, and protecting all kinds of your flexible endoscopes.

This product is certified as a medical device in the European Union under the Medical Device Directive 93/42/EEC by SGS CE0120, exclusively for the intended use of disinfects the flexible endoscope. Other non medical uses ascribed to this device are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those.

Multi-Channel jetting system rotated in a 360° arc



Automatic Water Line Sterilization (self disinfection) cycle



UV lamp disinfection system in the water tank



SPECIFICATION

Leakage Tester (water Type) / Operating Procedure

| | | |
|---------------------|--|----------------------------------|
| Time | Washing | 10~50sec, 1min~5min |
| | Disinfection | 1 ~ 9 min, 10 ~ 50 min, 1 ~ 9 hr |
| | Enzimatik | 10 sec |
| | Alcohol | 5 sec |
| Program (Automatic) | A mode : Washing-Detergent-Washing-Disinfectant-Washing-Alcohol-Air B mode : Washing-Detergent-Washing-Air C mode : Washing-Disinfectant-Washing-Air D mode : same as A mode (disinfectant is discharged after 1 cycle) | |
| TANK CAPACITY | Water | 20 Liters |
| | Disinfectant | 20 Liters |
| | Detergent : | 1 Liters |
| | Alcohol : | 1 Liters |
| COUPLERS | 6 Outlets in cleaning Base (2 Scopes) | |
| Discharge | 1 Outlet (Compulsory discharging on rear panel) | |
| | 3 Outlets (in the front door) | |
| Power | AC 220/230V, 50/60Hz | |
| Weight | 640 × 620 × 1050mm (W × D × H) | |
| Dimension | 86 kg. | |

FEATURES

- Simultaneous dual-scope cleaning and disinfection
- Leakage tester
- Perfect washing and disinfection by 3 powerful motors
- Multi-channel jetting system rotated in a 360 arc
- Advanced water filtration system
- UV lamp water disinfection system
- Thermal control system for disinfection
- Self disinfection system
- Alcohol purging system
- Dynamic printing system
- Deodorization system
- Automated lid opening system