




# Disposable Polypectomy Snare

Everything in a safe grip.





Removal of polyps or other excess tissues in the digestive tract used with diathermic energy.



## Special features

- Made by high-strength medical braided wire, offering a precise & quick cutting properties, as well as an outstanding electric coagulation.
- I Loop rotates synchronously by rotating the 3-ring handle, highly increase efficiency.
- I Standard power connection, compatible with all main high-frequency devices on the market.
- I Ergonomic design of 3-ring handle, easier to hold and use.

|   | Model                  | Diameter of the outer tube (mm) | Length (mm) | Trap open width (mm) |
|---|------------------------|---------------------------------|-------------|----------------------|
|  | <b>The Oval</b>        |                                 |             |                      |
|   | FMD-QTQ-TY-18x500x10   | 1.8                             | 500         | 10                   |
|   | FMD-QTQ-TY-18x1200x15  | 1.8                             | 1200        | 15                   |
|   | FMD-QTQ-TY-23x1600x20  | 2.3                             | 1600        | 20                   |
|   | FMD-QTQ-TY-23x2300x25  | 2.3                             | 2300        | 25                   |
|  | <b>Hexagon -shape</b>  |                                 |             |                      |
|   | FMD-QTQ-LJ-18x500x10   | 1.8                             | 500         | 10                   |
|   | FMD-QTQ-LJ--18x1200x15 | 1.8                             | 1200        | 15                   |
|   | FMD-QTQ-LJ--23x1600x20 | 2.3                             | 1600        | 20                   |
|   | FMD-QTQ-LJ--23x2300x25 | 2.3                             | 2300        | 25                   |