

double ends cleaning brush



## Disposable Cleaning Brush

### Intended use


- Cleaning of endoscope working and vapor channels




### vapor cleaning brush

- Used for cleaning vapor channel

## Cleaning Brush



Three-brushed cleaning brush: Disposable; equipped with three brushes which largely improve the cleaning efficiency and cleanliness



Double-brushed cleaning brush: Disposable.

Model	O.D.(mm)	Length(mm)	Endoscope Channel (mm)	Characteristics
JHY-CB-06-180-D	6	1800	2.8/3.2	With innovative three brushes, it has better cleaning result; Disposable ; Nylon brush ; Plastic shaft ;Three brushes;
JHY-CB-06-230-D	6	2300	2.8/3.2	
JHY-CB-03-120-D	3	1200	2.0	
JHY-CB-03-180-D	3	1800	2.0	Disposable ;Nylon brush ; Plastic shaft ;Double brushes;
JHY-CB-06-180-C	6	1800	2.8/3.2	
JHY-CB-06-230-C	6	2300	2.8/3.2	
JHY-CB-03-120-C	3	1200	2.0	Reusable ;Nylon brush ; Stainless steel shaft ; Single brush;
JHY-CB-03-180-C	3	1800	2.0	
JHY-CB-06-230-A	6	2300	2.8/3.2	
JHY-CB-03-180-A	3	1800	2.0	Reusable ;Nylon brush ; Stainless steel shaft ;Double brushes;
JHY-CB-06-230-B	6	2300	2.8/3.2	
JHY-CB-E	3	1800	2.0	Reusable ;Nylon brush ;Plastic handle;
JHY-CB-03-180-B	11/5	160	Vapor channel	