



Tek Kullanımlık Yabancı cisim Forceps

ENDOcept


Tanımlayıcı Özellikler



Model	O.D. (mm)	Boy (mm)	Endoskopi Kanalı (mm)	Karakteristik Özellikler
FMD-FG-23-160-B1	2.3	1600	≥2.8	Üç ayaklı, yuvarlak tırnaklı
FMD-FG-23-180-B1	2.3	1800	≥2.8	
FMD-FG-23-230-B1	2.3	2300	≥2.8	
FMD-FG-P-18 3C180	1.8	1800	≥2.0	
FMD-FG-P-24 3C230	2.4	2300	≥2.8	
FMD-FG-P-18 3R160	1.8	1600	≥2.0	
FMD-FG-P-24 3R230	2.4	2300	≥2.8	
FMD-FG-P-18 3R180	1.8	1800	≥2.0	
FMD-FG-P-24 3C180	2.4	1800	≥2.8	
FMD-FG-P-18 3C160	1.8	1600	≥2.0	
FMD-FG-P-18 3C120	1.8	1200	≥2.0	
FMD-FG-P-24 3R180	2.4	1800	≥2.8	
FMD-FG-P-18 3R120	1.8	1200	≥2.0	



Model	O.D. (mm)	Boy (mm)	Endoskopi Kanalı (mm)	Karakteristik Özellikler
FMD-FG-23-160-B2	2.3	1600	≥2.8	Dört Ayaklı, Yuvarlak Tırnaklı
FMD-FG-23-180-B2	2.3	1800	≥2.8	
FMD-FG-23-230-B2	2.3	2300	≥2.8	
FMD-FG-P-24 4C230	2.4	2300	≥2.8	
FMD-FG-P-24 4R180	2.4	1800	≥2.8	
FMD-FG-P-24 4C180	2.4	1800	≥2.8	
FMD-FG-P-24 4R230	2.4	2300	≥2.8	



Model	O.D. (mm)	Boy (mm)	Endoskopi Kanalı (mm)	Karakteristik Özellikler
FMD-FG-23-160-B3	2.3	1600	≥2.8	Beş Ayaklı, Yuvarlak Tırnaklı
FMD-FG-23-180-B3	2.3	1800	≥2.8	
FMD-FG-23-230-B3	2.3	2300	≥2.8	