



Disposable Biopsy Forceps

Quality and ergonomics perfectly combined.

The disposable biopsy forceps from stand out due to their high quality and ergonomic design. These features enable an easy and safe removal of tissue samples or polyps. Thanks to the first class manufacture, the biopsy forceps can be inserted very smoothly into the working channel and can be controlled precisely even with sharply angled endoscopes. The extra sharp jaws can be closed accurately and thus ensure a clean cut with minimal effort.

The extensive range offers the perfect solution for each application: whether with spike or without, with smooth or serrated jaws, with coating or without.

Ergonomic handle for optimum control.

For precise control of the jaw each biopsy forceps possess an ergonomically designed handle. It is particularly comfortable to hold, what makes the work easy and relaxed. Every the smallest rotation or movement of the handle is passed exactly to the jaw part so that this is controlled absolutely reliable.

Special features

- High-quality materials and manufacturing
- Simple and precise to use
- Sharp cutting edge for diagnostically conclusive biopsies
- Full closure of the cutting edges
- Special scissor design conserves the working channel
- Extensive product range

	Model	Diameter of jaws (mm)	Length (mm)
Uncoated without spike			
	FMD-QYQ-A-18X1050	1.8	1,050
	FMD-QYQ-A-18X1600	1.8	1,200
	FMD-QYQ-A-23X1600	1.8	1,600
	FMD-QYQ-A-23X1800	2.3	1,600
	FMD-QYQ-A-23X2300	2.3	1,800
		2.3	2,300
Uncoated with spike			
	FMD-QYQ-B-23X1600	2.3	1,600
	FMD-QYQ-B-23X1800	2.3	1,800
	FMD-QYQ-B-23X2300	2.3	2,300
Coated without spike			
	FMD-QYQ-C-18X1050	1.8	1,050
	FMD-QYQ-C-18X1200	1.8	1,200
	FMD-QYQ-C-18X1600	1.8	1,600
	FMD-QYQ-C-23X1600	2.3	1,600
	FMD-QYQ-C-23X1800	2.3	1,800
	FMD-QYQ-C-23X2300	2.3	2,300
Coated with spike			
	FMD-QYQ-D-23X1600	2.3	1,600
	FMD-QYQ-D-23X1800	2.3	1,800
	FMD-QYQ-D-23X2300	2.3	2,300
With alligator uncoated without spike			
	FMD-QYQ-E-23X1600	2.3	1,600
	FMD-QYQ-E-23X1800	2.3	1,800
	FMD-QYQ-E-23X2300	2.3	2,300
With alligator uncoated with spike			
	FMD-QYQ-F-23X1600	2.3	1,600
	FMD-QYQ-F-23X1800	2.3	1,800
	FMD-QYQ-F-23X2300	2.3	2,300
With alligator coated without spike			
	FMD-QYQ-G-23X1600	2.3	1,600
	FMD-QYQ-G-23X1800	2.3	1,800
	FMD-QYQ-G-23X2300	2.3	2,300
With alligator coated with spike			
	FMD-QYQ-B-23X1600	2.3	1,600
	FMD-QYQ-B-23X1800	2.3	1,800
	FMD-QYQ-B-23X2300	2.3	2,300