




Disposable Grasping Forceps C With

Everything in a safe grip.

Usage : It is used to grasp and remove foreign body from the alimentary tract in human body in combination with endoscope.

Special features

- Integral rotation handle
- With god flexibility and net design, it can accurately and firmiy grasp the forgign body.
- Singleuse, avoiding cross conlaminaior
- ESDand EMR's indispensable pariner.ihe clots,polyps, tissues can be thoroughiy retrisved.



The design of the handle is rotatable , the net shap grasper design avoid the lost ol samples, deeper elastic nylon net , which will efficienitly filter out the natural tract fluid,it is good for grasp and control .

Model	Diameter of the outhur tube (mm)	Length (mm)	Actual Opening width
FMD-YWQ-TY-18x1600x20	1.8	1600	15
FMD-YWQ-TY-23x1600x20	2,3	1600	15
FMD-YWQ-TY-23x1600x30	2.3	1600	25
FMD-YWQ-TY-23x1600x40	2.3	1600	35
FMD-YWQ-TY-23x1600x20	2,3	2300	15
FMD-YWQ-TY-23x1600x30	2.3	2300	25
FMD-YWQ-TY-23x1600x40	2.3	2300	35

