


Tek Kullanımlık Yabancı Cisim Forsepisi (Fileli)

ENDOcept

Kullanım Amacı

Sindirim kanalından yabancı cisimlerin yakalanmasını sağlar.



Tanımlayıcı Özellikler

Model	O.D.(mm)	Boy (mm)	Endoskopi Kanalı (mm)	Ağ büyüklüğü (mm)
FMD-FG-18-120-C1	1.8	1200	≥2.0	20*50
FMD-FG-18-120-C2	1.8	1200	≥2.0	25*55
FMD-FG-18-120-C3	1.8	1200	≥2.0	30*60
FMD-FG-18-160-C1	1.8	1600	≥2.0	20*50
FMD-FG-18-160-C2	1.8	1600	≥2.0	25*55
FMD-FG-18-160-C3	1.8	1600	≥2.0	30*60
FMD-FG-25-230-C1	2.5	2300	≥2.8	20*50
FMD-FG-25-230-C2	2.5	2300	≥2.8	25*55
FMD-FG-25-230-C3	2.5	2300	≥2.8	30*60
FMD-FG-25-230-C4	2.5	2300	≥2.8	40*75
FMD-FG-25-230-C5	2.5	2300	≥2.8	40*65
FMD-FG-W-24-352300	2.4	2300	≥2.8	35
FMD-FG-W-24-252300	2.4	2300	≥2.8	25

