



FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.

INFUSION SET SERIES



➤ Disposable Infusion Set with Latex Connector

Item Code: CMIF-1

Item Code: CMIF-1

Type: Luer slip or luer lock

Material: Medical grade PVC

Specification: 1. Drip chamber with air vent

2. Medical grade PVC tube, 120~180mm

3. Rubber drug-adding port

4. Polybag or blister packing

5. With hypodermic needle or with scalp vein set

6. Drip chamber: 48mm or 54mm, 20 drips / 1ml for adult; 60 drips / 1ml for pediatric

7. Gravity use only

INFUSION SET SERIES



➤ Disposable Infusion Set with Y-connector

Item Code: CMIF-2

Type: Luer slip or luer lock

Material: Medical grade PVC

Specification:


1. Drip chamber with air vent
2. Medical grade PVC tube, 120~180mm
3. Y drug-adding port
4. Polybag or blister packing
5. With hypodermic needle or with scalp vein set
6. Drip chamber: 48mm or 54mm, 20 drips / 1ml for adult; 60 drips / 1ml for pediatric
7. Gravity use only





FAMED SAĞLIK TIBBİ CİHAZLAR
İTH. İHRC. SAN. ve TİC. LTD. ŞTİ.

BURETTE INFUSION SET



➤ Burette Infusion Set

Item Code: CMIFB

Size: 100ml, 150ml

60 drops = 1 ± 0.1 ml

Material: medical grade PVC or PETG

BURETTE INFUSION SET



> Blood Transfusion Set

Item Code: CMBT

Product details: 1. Pike with or without air vent
2. ABS or PP roller clamp
3. Latex bulb or Y connector
4. Luer slip or luer lock connector
5. Approximately 20 drops/ml.

With 200 micron filter and 18 G hypodermic needle



> Infusion Set for Pump

Item Code: CMIP

Material: medical grade PVC

Specification: 1. Small / big chamber with air vent
2. Medical grade PVC tube, 120~180mm
3. Polybag or blister packing
4. Drip chamber: 48mm or 54mm, 20 drips / 1ml
5. It can be used with the pump