



Disposable Cytology Brush

Innovational brush design, with no risk of drop-off
 Straight-shaped brush: easy to enter the depths of respiratory and digestive tract with
 U-shaped brush: obtain tissues from all directions with a large quantity

Endoscopy Cytology Brush

- **Usage:** Be used for cell tissue sampling from halitus tract and alimentary tract in the observation of endoscope.
- Brush of integration, without possibility of breaking off.
- With high smoothness, the outer tube is not easy to get bended
- With suitable stiffness of bristles, it is easy for cell tissue sampling
- Sterilized, only for disposable use.



Model	O.D.(mm)	Length(mm)	Endoscope Channel (mm)	Characteristics
FMD-BC-18-105-A	1.8	1050	φ2.0	Straight Shaped Brush
FMD-BC-18-120-A	1.8	1200	φ2.0	Straight Shaped Brush
FMD-BC-23-160-A	2.3	1600	φ2.8	Straight Shaped Brush
FMD-BC-23-180-A	2.3	1800	φ2.8	Straight Shaped Brush
FMD-BC-23-230-A	2.3	2300	φ2.8	Straight Shaped Brush
FMD-BC-18-105-B	1.8	1050	φ2.0	"U" Shaped Brush
FMD-BC-18-120-B	1.8	1200	φ2.0	"U" Shaped Brush
FMD-BC-23-160-B	2.3	1600	φ2.8	"U" Shaped Brush
FMD-BC-23-180-B	2.3	1800	φ2.8	"U" Shaped Brush
FMD-BC-23-230-B	2.3	2300	φ2.8	"U" Shaped Brush