

Table Top Ultrasonic Cleaning Units



4UT Pro D



TECHNICAL SPECIFICATIONS

Chassis	AISI 304 stainless steel	is made of certified stainless steel in accordance with health codex.
External Dimensions (MM)	300 x 255 x 230	
Tank Dimensions (MM)	240 x 135 x 150	
Water Tank Capacity (lt)	4	Stainless
Ultrasonic frequency (kHz)	40	
Ultrasonic Power (W)	60	
Transducer	1 pcs /2 piezoceramic	ø38MM 38,000 hours approx. Working life
Heater power (W)	200	Special production, protected against waterless operation.
Total maximum power (W)	260	
Timer (min)	0-59	Digital programmable, double digit display and buzzer.
Thermostat (° C)	0-75	Digital programmable, double digit display and buzzer.
Water drain	No	
Stainless Basket	Yes	

Usage Areas

- Hospital equipment
- Surgical instruments
- Veterinarians
- Dentist
- Disinfection applications
- Laboratories
- Optical sector
- Pharmaceutical industry

