

Table Top Ultrasonic Cleaning Units



12UT Pro D



TECHNICAL SPECIFICATIONS

Chassis	AISI 304 Stainless Steel	is made of certified stainless steel in accordance with health codex.
External Dimensions (MM)	315 x 290 x 320	
Tank Dimensions (MM)	300 x 240 x 200	
Water Tank Capacity (lt)	12	Stainless
Ultrasonic frequency (kHz)	40	
Ultrasonic Power (W)	180	
Transducer	3 pcs / 6 piezoceramic	ø38MM 38,000 hours approx. Working life
Heater power (W)	600	Special production, protected against waterless operation.
Total maximum power (W)	~780	
Timer (min)	0-59	Digital programmable, double digit display and buzzer.
Thermostat (° C)	0-75	Digital programmable, double digit display and buzzer.
Water drain	Yes	
Stainless Basket	Yes	

Usage Areas

- Hospital equipment
- Surgical instruments
- Veterinarians
- Dentist
- Disinfection applications
- Laboratories
- Optical sector
- Pharmaceutical industry

