

## 18 UT Pro D



### TECHNICAL SPECIFICATIONS

Chassis	AISI 304 Stainless Steel	is made of certified stainless steel in accordance with health codex.
External Dimensions (MM)	330 x 350 x 500	
Tank Dimensions (MM)	315 x 320 x 200	
Water Tank Capacity (lt)	14	Stainless
Ultrasonic frequency (kHz)	28	
Ultrasonic Power (W)	240	
Transducer	4 pcs / 8 piezoceramic	ø38MM 38,000 hours approx. Working life
Heater power (W)	800	Special production, protected against waterless operation.
Total maximum power (W)	~1040	
Timer (min)	0-59	Digital programmable, double digit display and buzzer.
Thermostat (° C)	0-75	Digital programmable, double digit display and buzzer.
Water drain	Yes	
Stainless Basket	Yes	

### Usage Areas

- Hospital equipment
- Surgical instruments
- Veterinarians
- Dentist
- Disinfection applications
- Laboratories
- Optical sector
- Pharmaceutical industry