

Table Top Ultrasonic Cleaning Units



28 UT Pro D

TECHNICAL SPECIFICATIONS

Chassis	AISI 304 Stainless Steel	is made of certified stainless steel in accordance with health codex.
External Dimensions (MM)	500 x 300 x 200	
Tank Dimensions (MM)	515 x 315 x500	
Water Tank Capacity (lt) / Water capacity (lt)	28/20	Stainless
Ultrasonic frequency (kHz)	28	
Ultrasonic Power (W)	360	
Transducer	6 pcs / 12 piezoceramic	ø38MM 38,000 hours approx. Working life
Heater power (W)	1000	Special production, protected against waterless operation.
Total maximum power (W)	~1360	
Timer (min)	0-59	Digital programmable, double digit display and buzzer.
Thermostat (° C)	25-75	Digital programmable, double digit display and buzzer.
Water drain	Yes	
Stainless Basket	Yes	

Usage Areas

- Hospital equipment
- Surgical instruments
- Veterinarians
- Dentist
- Disinfection applications
- Laboratories
- Optical sector
- Pharmaceutical industry

