



BIORAD®

THERMO CLEAN %21

Hemodialysis Device
Disinfection



FORMS OF PACKAGE

| Description | Package |
|-------------|---------|
| 5 L PE drum | 4x5 L |

CHEMICAL COMPOSITION:

Citric Acid, Monohydrate, Malic Acid, Lactic Acid and auxiliary agents

PRODUCT SPECIFICATIONS:

THERMO CLEAN %21 provides cleaning, decalcification and disinfection of organic and inorganic waste that forms on the generator and lines of haemodialysis equipment during haemodialysis process. It is compliant with stainless steel, plastic materials and ceramic surfaces. It should not contact with plexiglass and aluminium surfaces.

INTENDED USE:

Thermo Clean %21 is compliant with stainless steel, glass, plastic materials and ceramic surfaces. Do not use it on sensitive surfaces such as plexiglass and aluminium that may react with acids.

USAGE:

Thermo Clean %21 is used by automatic vacuum through disinfection pipette of the haemodialysis achines. It is recommended to use it at rate of 1/120. pH=1-3

ANTIMICROBIAL CHARACTERISTICS:

Bactericidal, Fungicidal, Virucidal, Micobactericidal.Tuberculocidal