

# DSC ENDO 10

Safe, Robust Endoscope Storage Cabinet



FAMED SAĞLIK TIBBİ CİHAZLAR İTH.İHRC. ve SAN. TİC. LTD. ŞTİ.



## FEATURES

- Endoscopic cabinet is designed for drying and sterilizing the cleaned endoscope by fan heater and UV lamp.
- This is designed to keep high-priced endoscope safely and efficiently in endoscopy center and endoscope examination room and provide more convenience with user.
- Heater is equipped in interior to dry moisture after cleaning endoscope
- UV lamp is equipped to sterilize to upside and Timer to keep uniform temperature in chamber.
- All functions are operated automatically via sensor attached in door side when door and close.
- SIZE(LxDxH,mm):655\*655\*2,225
- UV lamp 20W(life time 5,00hr)
- Power: AC220V 50/60Hz
- Control box : On& Off power switch, Timer (30min.)

- 10 HOLDERS
- ROTATING RACK
- UV-C RADIATOR
- DRYING CABINET
- TIMER CONTROLLER
- CLEAR VIEW
- COMPACT DESIGN

