

DSC ENDO -8

Safe, Robust Endoscope Storage Cabinet





FAMED SAĞLIK TIBBİ CİHAZLAR İTH.İHRC. ve SAN. TİC. LTD. ŞTİ.



*Medical equipment manufacturer,
Endoscope cleaner, Sterilizer,
therapy laser, Disinfectant and more*

- ▶ 8 Holders
- ▶ UV-C Radiator
- ▶ Automatic operation by door sensor
- ▶ Sliding door
- ▶ Lockable cabinet
- ▶ SUS304 material
- ▶ Air circulation fan



FEATURES

- Endoscopic cabinet is designed for drying and sterilizing the cleaned endoscope by fan heater and UV lamp
- This is designed to keep high-priced endoscope safely and efficiently in endoscopy center and endoscope examination room and provide more convenience with user.