



BIORAD®

THERMO CLEAN %30

Hemodialysis Device
Disinfection



FORMS OF PACKAGE

Description	Package
5 L PE drum	4x5 L

CHEMICAL COMPOSITION:

%30 Citric Acid, Monohydrate and auxiliary agents

PRODUCT SPECIFICATIONS:

THERMO CLEAN %30 provides cleaning, decalcification and disinfection of organic and inorganic waste that forms on the generator and lines of haemodialysis equipment during haemodialysis process. It is compliant with stainless steel, plastic materials and ceramic surfaces. It should not contact with plexiglass and aluminium surfaces.

INTENDED USE:

Thermo Clean %30 is compliant with stainless steel, glass, plastic materials and ceramic surfaces. Do not use it on sensitive surfaces such as plexiglass and aluminium that may react with acids.

USAGE:

Thermo Clean %30 is used by automatic vacuum through disinfection pipette of the haemodialysis machines. Observe recommendations of the machine manufacturer. pH=1-3

ANTIMICROBIAL CHARACTERISTICS:

Bactericidal, Fungicidal, Virucidal, Micobactericidal, Tuberculocidal